

## Evok 800 - Key Features & Benefits

### Behind the Ear (BTE) style

This is the latest design BTE style hearing aid and is small and discreet. It adapts to your environment for crystal clear sound quality and filters out background noise to focus on speech in any setting. Powerful and completely tailored to you, this hearing aid will transform the way you hear.

### Digital / Programmable

Fully programmable hearing aid that is programmed to your hearing loss and preferences. It is tailored to you. Features 4 programs for different listening environments:

1. Everyday General Use
2. Background noise situations (restaurants etc)
3. Phone Use and talking in small group situations
4. Music and TV

### Adaptive Directional Microphones

Two microphones mean that when the relevant program is selected the aid detects the person speaking in front of the wearer and automatically configures the microphones to focus in the speakers direction, whilst reducing noise from other directions. This helps improve wearer comfort in noisy situations.

### Hydrocoat

This helps your hearing aids resist moisture, sweat, oil, and corrosion. Hydrocoat is an advanced micron-level sealant.

### Volume/Memory Retention

Even when turned off the hearing aid will remember your volume and memory settings next time you turn on.

### Adaptive Feedback Shield

This helps eliminate feedback and whistling. While all styles of aids have limits to the amount of feedback they can cancel, this is maximised so it can “seek and destroy” feedback across a wide band of frequencies. Feedback-Shield utilises an adaptive phase canceller, which means if your surroundings change and the signal causing feedback changes, then the aid changes.

### Tri-Mode Layered Noise Reduction

Emphasise speech and improve clarity in challenging situations (noisy environment) and helps make listening natural and comfortable.

### Dynamic Speech Enhancement (DSE)

Designed to enhance speech understanding in difficult environments and enable you to listen to both a soft and loud voice during the same conversation. Sudden changes to signal levels are automatically adjusted for by the micro-processor.

### Telecoil

**Built in Telecoil means you can connect to phones with Telecoil and anywhere there is a Hearing Loop system like movie theatres, churches, concert halls, universities, museums and airports etc. This allows your hearing aids to wirelessly stream the sound of a speaker's voice or the music directly into your ears – without any background noise.**

**Voice Prompts and Low Battery Indicator**

**The Evok 800 has voice prompts for program changes and low battery warning.**

